

Adjustable tables

# Solu-Tables

TAB-C, TAB-R & TAB-C Series

FOCUSED ON WHAT MATTERS



SPECIALTY  
HEALTHCARE  
SEATING



TAB4T MODEL  
CLOVER - 4 PLACES

ADJUSTABLE TABLES

# Solu-Tables

## Adjustable tables

The Solu-Tables are designed with ease of use in mind. Their sliding system makes adjusting the trays effortless and safe, providing great stability when sliding with a food tray on top. Thanks to the levelers on each leg, leveling the table is a breeze and leaves no marks on the floor. The base is made of rectangular tubes to reduce trip hazards and provide optimal clearance for the resident's feet.

The concave-shaped trays allow wheelchairs to get very close to the table, maximizing space for food service. The leg is designed to offer additional clearance between the floor and the base, making floor cleaning easier. A hidden wheel helps staff move the table with minimal effort and risk of injury.



# Resistant, stable and safe tables

For over 30 years, LPA Medical has been manufacturing medical equipment designed to help people experiencing autonomy loss. Entirely made in Quebec, our dining room tables help facilitate the work of nursing staff as well as the comfort of residents.

## TECHNICAL SPECIFICATIONS

	TAB1R	TAB2R	TAB3T TAB4T	TAB4C	TAB4R	TAB6R
Tray width	66 cm / 26"		64 cm / 25"	63,5 cm / 25"	66 cm / 26"	
Tray length	56 cm / 22"		64 cm / 25"	89 cm / 35"	56 cm / 22"	
Total width (overall)	66 cm / 24"	66 cm / 26"	132 cm / 52"	Diameter 131 cm / 52"	137 cm / 54"	249 cm / 98"
Total length (overall)	61 cm / 24"	117 cm / 46"			117 cm / 46"	117 cm / 46"

Please note that the technical specifications are subject to change without notice.



## Available tray finishes

Tray finishes are available in 4 standard colors designed to complement most existing furniture arrangements. Having these 4 standard finishes helps guarantee fast lead times. Should you need a specific color, a wide variety of options exists and we'll be happy to assist you with your selection. Contact us for more information.



MAHOGANY  
7008-58



WALNUT  
3485-58



NATURAL OAK  
7737-58



LIGHT MAPLE  
86992-58

# 7 models of adjustable tables



TAB1R MODEL  
RECTANGULAR – 1 PLACE



TAB2R MODEL  
RECTANGULAR – 2 PLACES



TAB3T MODEL  
CLOVER – 3 PLACES



TAB4T MODEL  
CLOVER – 4 PLACES



TAB4R MODEL  
RECTANGULAR – 4 PLACES



TAB6R MODEL  
RECTANGULAR – 6 PLACES



TAB4C MODEL  
ROUND – 4 PLACES

Visit our website or contact us for more details about our range of products.



[lpamedical.com](http://lpamedical.com)

2527 Dalton avenue, Quebec (Qc) G1P 3S6 CANADA  
T 418-681-1313 1-800-663-4863 F 418-681-4488

Edited on 2023-12