

patient care recliner

Harmony

C Series

FOCUSED ON WHAT MATTERS



SPECIALTY
HEALTHCARE
SEATING



HARMONY C-207

Shown in optional
Silvertex Plata vinyl



HARMONY C-207

Shown in optional
Silvertex Sapphire vinyl

COMFORT AND POSITIONING

Harmony

patient care recliner



Optional Silvertex vinyl containing silver nitrate*
Available in several colours, it is an ideal material
to help combat the spread of infections

* Silver nitrate is recognized for its antibacterial properties

Removable high resilience foam cushions
for easy cleaning and superior comfort

Two swivel wheels with brakes in
a diagonal configuration, a directional
wheel and a single swivel wheel standard

Ergonomic, Attractive and Sturdy

Seat tilting system allows
the occupant to push back
and recline easily

Upper backrest with push bar
(helps the caregiver
tilt the chair)



The ideal combination of comfort and function

Providing height adjustability, versatility and excellent value, the Harmony patient care recliner can be used for taking blood samples, specific treatments, general care and prolonged sitting. It is suitable for both short and long-term care establishments.



C-C 1



C-C 2



C-C 4




C-C6



C-C 8



C-C 9



C-C 10



C-C 11



RB-18

OPTIONS AND ACCESSORIES

ITEM	DESCRIPTION
C-C1	Left or right lateral tray
C-C2	Adjustable elbow supports
C-C3	Adjustable IV pole, rear installation
C-C3A	Adjustable IV pole, front installation
C-C3B	Adaptor for IV pole with 1/2" round tube
C-C4	Elevating footrest support
C-C5	Set of four (4) swivel wheels with brakes
C-C6	Self-storing, removable footrests
C-C7	Side bolster cushions
C-C8	Opaque activity tray
C-C8-2	Activity tray with Lexan window
C-C9	Lateral support bracket for storing activity tray
C-C10	PRE-CURE® PRO cushion allows the insertion of speciality pressure management gels/foams
C-C11	Inclined armrest for blood sample
AB-08H	Contoured headrest
AC-04	Leatherette padded tray cover
AE-02	Catheter drainage bag hook
RA-03	Lateral tray straps
RA-03A	Wrap-around tray strap
RA-19	Universal perineal abductor kit with black fixing mat
RA-23	4" (10 cm) lateral trunk supports with backrest canvas
RA-24	6" (15 cm) lateral trunk supports with backrest canvas
RB-18	Cervical support

Vinyl upholstery standard. Crypton fabric and Silvertex vinyl available as options. Meets CA TB117-2013 fire standards.

TECHNICAL SPECIFICATIONS

	C-207-20
Overall width	28 1/2" / 72 cm
Upright length	36 1/2" / 93 cm
Reclined length	70" / 178 cm
Upright height	46 1/2" / 118 cm
Seat/floor height	Adjustable from 17 1/2" to 20 1/2" (44 to 52 cm)
Seat width	20 1/4" / 51 cm
Seat depth	Adjustable from 17" to 18 1/2" (43 to 47 cm)
Seat/backrest angle	Adjustable from 95 to 108 degrees
Weight w/o accessories	80 lbs. / 36 kg
Maximum weight capacity	250 lbs. / 114 kg

Please note that technical specifications are subject to change without notice.

our complete product range

MAINTAINING PATIENT AUTONOMY

THERA-GLIDE®
WOOD



THERA-GLIDE®
METAL



THERA-GLIDE®
TRANSITIONAL



DYN-ERGO



NOMAD



VERSATILE



EVOLUTION



UTILITY



COMFORT AND POSITIONING

COMP-POSITION



HARMONY




OPTIMUM



ACCESSORIES

PRE-CURE® PRO
CUSHION



ANTI-SLIP
SURFACES



Visit our website or contact us for more details about our range of products.

lpa+

lpamedical.com

2527 Dalton avenue, Quebec (Qc) G1P 3S6 CANADA
T 418-681-1313 1-800-663-4863 F 418-681-4488