

mobility chairs

# Dyn-Ergo and Evolution S Series



FOCUSED ON WHAT MATTERS

SPECIALTY  
HEALTHCARE  
SEATING



DYN-ERGO S-150

Shown with 13" (SA-30)  
wheels and optional Silvertex  
Turquoise vinyl



EVOLUTION S-250

Shown in standard Naugasoft Cocoa vinyl  
with removable/foldable footrests (SA-15),  
optional cervical support (RB-18)  
and thoracic support (SA-20)

MAINTAINING PATIENT AUTONOMY

# Dyn-Ergo and Evolution

mobility chairs

## Standard adjustments

Optional Silvertex vinyl containing silver nitrate\*  
Available in several colours, it is an ideal material  
to help combat the spread of infections

\* Silver nitrate is recognized for its antibacterial properties

PRE-CURE® seat cushion composed of high resilience foam  
and visco gel memory foam  
This cushion is standard and is designed  
to help prevent pressure sores

Adjustable Backrest  
(95° to 115°)


3 possible  
armrest heights  
(8", 8½" and 9½"  
/ 20, 22 and 24 cm)

Seat depth  
(13½" to 18" / 35 to 46 cm)

Seat height  
S-150: 13½" to 18½" / 34 to 47 cm  
S-250: 16" to 21" / 41 to 53 cm

20" central wheels for  
instinctive mobility  
(Model S-250)

The maximum tilt  
of 34° and the cylinder's  
8 adjustable locking positions  
make it possible to adapt it  
to most of the user's needs



# Our mobility chairs...

- + Meet mobility and positioning needs and are an advantageous alternative to wheelchairs and/or bedroom armchairs.
- + Provide greater freedom of movement for patients at risk of falling and help reduce agitation for those with cognitive impairment.
- + Help maintain patient autonomy and simplify reassignment thanks to their multiple adjustments and accessories.
- + Facilitate social acceptance through their user-friendly appearance and comfort.



SA-10



SA-15




SA-30



SA-40



SA-20



SA-60



SA-130



SA-50



RA-23

## OPTIONS AND ACCESSORIES

ITEM	DESCRIPTION	Dyn-Ergo	Évolution
SA-10	Swing-away and removable footrests	0	—
SA-15	Swing-away and removable footrests	—	0
SA-20	Structural lateral supports	0	0
SA-30	13" (33 cm) rear wheels with brakes	0	—
SA-40	Activity tray	0	0
SA-40-2	Activity tray with Lexan window	0	0
SA-50	Cervical support	0	0
SA-60	PRE-CURE® PRO cushion enabling insertion of specialty pressure management gels/foams	0	0
SA-70	Heel support loops	0	0
SA-75	Height-adjustable IV pole	0	0
SA-80	Adjustable Seat depth 15 1/2" to 20 3/4" (39 to 53 cm)	0	0
SA-90	Seat height increase adjustment option from 16 1/2" to 21 1/2" (42 to 55 cm). Incompatible with option SA-30	0	0
SA-120	Full-width padded footrest (requires option SA-10 or SA-15)	0	0
SA-130	Adjustable elevating leg rest (sold individually or as a pair)	0	0
SA-215	Adjustable tension backrest	0	0
RA-03	Lateral tray straps	0	0
RA-03A	Wrap-around tray strap	0	0
RA-19	Perineal abductor with universal black fixing mat	0	0
RA-23	4" (10 cm) adjustable lateral trunk supports	0	0
RA-24	6" (15 cm) adjustable lateral trunk supports	0	0
AB-08-H	Contoured headrest	0	0
AE-02	Catheter drainage bag hook	0	0

+: Standard 0: Optional —: Not available

Vinyl upholstery standard. Crypton fabric and Silvertex vinyl available as options. Meets CA TB117-2013 fire standards.

## TECHNICAL SPECIFICATIONS

	Dyn-Ergo (S-150)					Evolution (S-250)				
	16	18	20	22	24	16	18	20	22	24
Maximum width	58" / 147 cm					57" / 145 cm				
Minimum width	31" / 79 cm					32" / 81 cm				
Overall width	23" / 58 cm	25 1/2" / 65 cm	27" / 69 cm	29" / 74 cm	30 1/2" / 77 cm	25 1/2" / 65 cm	28" / 71 cm	29 1/2" / 75 cm	31 1/2" / 80 cm	33" / 84 cm
Overall width with option SA-30	27" / 69 cm	29 1/2" / 75 cm	31" / 79 cm	33" / 84 cm	34 1/2" / 88 cm	—	—	—	—	—
Maximum height	54" / 137 cm					54 1/2" / 138 cm				
Minimum height	36 1/2" / 93 cm					38" / 97 cm				
Seat width	15 1/2" / 39 cm	18" / 46 cm	19 1/2" / 50 cm	21 1/2" / 55 cm	23" / 58 cm	15 1/2" / 39 cm	18" / 46 cm	19 1/2" / 50 cm	21 1/2" / 55 cm	23" / 58 cm
Seat depth	Adjustable from 13 3/4" to 18" (35 to 46 cm)					Adjustable from 16" to 21" / 41 to 53 cm				
Seat height	Adjustable from 13 1/2" to 18 1/2" / 34 to 47 cm					Adjustable from 16" to 21" / 41 to 53 cm				
Armrest height	Adjustable from 8" to 9 1/2" / 20 to 24 cm					Adjustable from 95 to 115 degrees				
Backrest tilt	Adjustable from 95 to 115 degrees					Adjustable from 95 to 115 degrees				
Maximum seat tilt	34 degrés					23 degrés				
Weight without accessories	46 lbs. / 21 kg	49 lbs. / 22 kg	52 lbs. / 24 kg	62 lbs. / 28 kg	65 lbs. / 29 kg	70 lbs. / 32 kg	73 lbs. / 33 kg	76 lbs. / 34 kg	86 lbs. / 39 kg	89 lbs. / 40 kg
Maximum weight capacity	250 lbs. / 113 kg			400 lbs. / 181 kg		250 lbs. / 113 kg			400 lbs. / 181 kg	



# our complete product range

## MAINTAINING PATIENT AUTONOMY



THERA-GLIDE®  
WOOD



THERA-GLIDE®  
METAL



THERA-GLIDE®  
TRANSITIONAL



THERA-GLIDE®  
OASIS



DYN-ERGO



EVOLUTION



NOMAD



VERSATILE



UTILITY

## COMFORT AND POSITIONING



COMP-POSITION



HARMONY



OPTIMUM

## ACCESSORIES



PRE-CURE® PRO  
CUSHION



ANTI-SLIP  
SURFACES

Visit our website or contact us for more details about our range of products.

lpa+

[lpamedical.com](http://lpamedical.com)

2527 Dalton avenue, Quebec (Qc) G1P 3S6 CANADA  
T 418-681-1313 1-800-663-4863 F 418-681-4488