

positioning tilt-in-space chair

# Comp-Position

0 Series

FOCUSED ON WHAT MATTERS

lpa+

SPECIALTY  
HEALTHCARE  
SEATING



COMP-POSITION 0-7

Shown in Silvertex Ice vinyl  
with removable armrests (0A-03)  
and 13" optional wheels (0E-01)



COMP-POSITION 0-07

Shown with standard  
Naugasoft Sapphire vinyl  
and optional removable  
armrests (0A-03)

COMFORT AND POSITIONING

# Comp-Position

positioning tilt-in-space chair




Standard articulated headrest  
for better head and neck support  
and alignment

PRE-CURE® seat cushion composed  
of high resilience foam and visco gel  
memory foam

This cushion is standard and is designed  
to help prevent pressure ulcers


Various adjustments:

- Seat width
- Floor-to-seat height
- Backrest width



2 simple functions

- Tilt-in-space
- Leg rests with automatic  
length adjustment



13" optional  
rear wheels (code 0E-01)  
with central brake




# Adjustable and easily reallocated

- + Meets the highest standards of positioning, comfort and security.
- + Suitable for different morphologies and posture problems.
- + Enables safe transfers both independently and with the help of a patient lift.
- + Makes it easier for caregivers to reduce the risk of accidents at work.
- + Reduces inventory and costs as its adjustments enable easy reassignment.



OA-03



OA-04



OB-02



OB-05



OB-07



OC-01



OC-04



OD-02



OE-01

## OPTIONS AND ACCESSORIES

| ITEM   | DESCRIPTION                                                   |
|--------|---------------------------------------------------------------|
| AE-02  | Catheter drainage bag hook                                    |
| OA-01  | 1 1/2" side bolsters                                          |
| OA-03  | Removable armrests                                            |
| OA-04  | PRE-CURE® PRO cushion                                         |
| OA-06  | Additional seat cushion                                       |
| OB-02  | Adjustable tension backrest                                   |
| OB-03  | 4" [10 cm] lateral trunk support                              |
| OB-04  | 6" [15 cm] lateral trunk support                              |
| OB-05  | Adjustable structural lateral supports                        |
| OB-06  | Adjustable shoulder supports                                  |
| OB-07  | Gas cylinder assisted back angle adjustment 98 to 135 degrees |
| OB-08  | Additional backrest cushion                                   |
| OC-01  | Opaque activity tray                                          |
| OC-02  | Activity tray with Lexan window                               |
| OC-03  | Padded leatherette tray cover                                 |
| OC-04  | Hemiplegic armrest pad                                        |
| OD-01  | Padded leatherette footboard cover                            |
| OD-02  | Lateral leg protectors                                        |
| OD-03  | Elevating footrest cushion                                    |
| OE-01  | 13" wheels with central brake                                 |
| OE-03  | IV pole                                                       |
| OE-012 | 13" wheel side brakes                                         |
| RA-03  | Lateral tray straps                                           |
| RA-03A | Wrap-around tray strap                                        |
| RA-19  | Universal perineal abductor kit with black fixing mat         |

Vinyl upholstery standard. Crypton fabric and Silvertex vinyl available as an option.  
Meets CA TB117-2013 fire standards..

## TECHNICAL SPECIFICATIONS

|                                                     | Série 0                                 |
|-----------------------------------------------------|-----------------------------------------|
| Overall width                                       | 25 3/4" / 65.4 cm                       |
| Overall depth (chair upright)                       | 34 1/2" / 87.6 cm                       |
| Overall depth (chair tilted, reclined and extended) | 65" / 165.1 cm                          |
| Overall height                                      | 53 1/4" / 135.3 cm                      |
| Floor-to-seat height                                | 19" to 21" / 48.3 cm to 53.3 cm         |
| Seat width between armrests                         | 16 1/4" to 19 1/4" / 41.3 cm to 48.9 cm |
| Seat depth                                          | 15 1/2" to 19 1/2" / 39.4 cm to 49.5 cm |
| Backrest height                                     | 22 1/2" / 57.2 cm                       |
| Backrest height with headrest                       | 34" / 86.4 cm                           |
| Length of leg rest (closed)                         | 16 1/2" to 18 1/2" / 41.9 cm to 47 cm   |
| Leg rest extension                                  | 4 3/4" to 5 3/4" / 12 cm to 14.6 cm     |
| Armrest height                                      | 7 1/4" to 12 1/4" / 18.4 cm to 31.1 cm  |
| Maximum tilt                                        | 45 degrees                              |
| Seat / backrest angle                               | 90 to 115 degrees                       |
| Total weight w/o accessories                        | 104 lbs. / 47 kg                        |
| Maximum weight                                      | 250 lbs. / 113.4 kg                     |

Please note that technical specifications are subject to change without notice.

# our complete product range

## MAINTAINING PATIENT AUTONOMY



THERA-GLIDE®  
WOOD



THERA-GLIDE®  
METAL



THERA-GLIDE®  
TRANSITIONAL



THERA-GLIDE®  
OASIS



DYN-ERGO



EVOLUTION



NOMAD



VERSATILE



UTILITY

## COMFORT AND POSITIONING



COMP-POSITION



HARMONY



OPTIMUM

## ACCESSORIES



PRE-CURE® PRO  
CUSHION



ANTI-SLIP  
SURFACES

Visit our website or contact us for more details about our range of products.

lpa+

[lpamedical.com](http://lpamedical.com)

2527 Dalton avenue, Quebec (Qc) G1P 3S6 CANADA  
T 418-681-1313 1-800-663-4863 F 418-681-4488