

positioning tilt-in-space chair

Comp-Position

O Series

FOCUSED ON WHAT MATTERS



SPECIALTY
HEALTHCARE
SEATING



COMP-POSITION O-7

Shown in Silvertex Ice vinyl
with removable armrests (OA-03)
and 13" optional wheels (OE-01)



COMP-POSITION O-07

Shown with standard
Naugasoft Sapphire vinyl
and optional removable
armrests (OA-03)

COMFORT AND POSITIONING

Comp-Position

positioning tilt-in-space chair




Standard articulated headrest for better head and neck support and alignment

PRE-CURE® seat cushion composed of high resilience foam and visco gel memory foam

This cushion is standard and is designed to help prevent pressure ulcers


Various adjustments:

- Seat width
- Floor-to-seat height
- Backrest width



2 simple functions

- Tilt-in-space
- Leg rests with automatic length adjustment




13" optional rear wheels (code 0E-01) with central brake

Adjustable and easily reallocated

- + Meets the highest standards of positioning, comfort and security.
- + Suitable for different morphologies and posture problems.
- + Enables safe transfers both independently and with the help of a patient lift.
- + Makes it easier for caregivers to reduce the risk of accidents at work.
- + Reduces inventory and costs as its adjustments enable easy reassignment.



OA-03



OA-04



OB-02



OB-05



OB-07



OC-01



OC-04



OD-02



OE-01

OPTIONS AND ACCESSORIES

ITEM	DESCRIPTION
AE-02	Catheter drainage bag hook
OA-01	1 1/2" side bolsters
OA-03	Removable armrests
OA-04	PRE-CURE® PRO cushion
OA-06	Additional seat cushion
OB-02	Adjustable tension backrest
OB-03	4" [10 cm] lateral trunk support
OB-04	6" [15 cm] lateral trunk support
OB-05	Adjustable structural lateral supports
OB-06	Adjustable shoulder supports
OB-07	Gas cylinder assisted back angle adjustment 98 to 135 degrees
OB-08	Additional backrest cushion
OC-01	Opaque activity tray
OC-02	Activity tray with Lexan window
OC-03	Padded leatherette tray cover
OC-04	Hemiplegic armrest pad
OD-01	Padded leatherette footboard cover
OD-02	Lateral leg protectors
OD-03	Elevating footrest cushion
OE-01	13" wheels with central brake
OE-03	IV pole
OE-012	13" wheel side brakes
RA-03	Lateral tray straps
RA-03A	Wrap-around tray strap
RA-19	Universal perineal abductor kit with black fixing mat

Vinyl upholstery standard. Crypton fabric and Silvertex vinyl available as an option.
Meets CA TB117-2013 fire standards..

TECHNICAL SPECIFICATIONS

	Série 0
Overall width	25 3/4" / 65.4 cm
Overall depth (chair upright)	34 1/2" / 87.6 cm
Overall depth (chair tilted, reclined and extended)	65" / 165.1 cm
Overall height	53 1/4" / 135.3 cm
Floor-to-seat height	19" to 21" / 48.3 cm to 53.3 cm
Seat width between armrests	16 1/4" to 19 1/4" / 41.3 cm to 48.9 cm
Seat depth	15 1/2" to 19 1/2" / 39.4 cm to 49.5 cm
Backrest height	22 1/2" / 57.2 cm
Backrest height with headrest	34" / 86.4 cm
Length of leg rest (closed)	16 1/2" to 18 1/2" / 41.9 cm to 47 cm
Leg rest extension	4 3/4" to 5 3/4" / 12 cm to 14.6 cm
Armrest height	7 1/4" to 12 1/4" / 18.4 cm to 31.1 cm
Maximum tilt	45 degrees
Seat / backrest angle	90 to 115 degrees
Total weight w/o accessories	104 lbs. / 47 kg
Maximum weight	250 lbs. / 113.4 kg

Please note that technical specifications are subject to change without notice.

our complete product range

MAINTAINING PATIENT AUTONOMY



THERA-GLIDE®
WOOD



THERA-GLIDE®
METAL



THERA-GLIDE®
TRANSITIONAL



THERA-GLIDE®
OASIS ROCKING



THERA-GLIDE®
OASIS STATIONNARY



DYN-ERGO



EVOLUTION



NOMAD



VERSATILE



UTILITY

COMFORT AND POSITIONING



COMP-POSITION



HARMONY



OPTIMUM

ACCESSORIES



PRE-CURE® PRO
CUSHION



ANTI-SLIP
SURFACES

Visit our website or contact us for more details about our range of products.



lpamedical.com

2527 Dalton avenue, Quebec (Qc) G1P 3S6 CANADA
T 418-681-1313 1-800-663-4863 F 418-681-4488

Edited on 2023-01