

patient care recliner

Harmony

C Series

FOCUSED ON WHAT MATTERS



SPECIALTY
HEALTHCARE
SEATING



HARMONY C-207

Shown in optional
Silvertex Plata vinyl



HARMONY C-207

Shown in optional
Silvertex Sapphire vinyl

COMFORT AND POSITIONING

Harmony

patient care recliner



Optional Silvertex vinyl
Available in several colours

Removable high resilience foam cushions
for easy cleaning and superior comfort

Three swivel wheels with brakes and
a front right directional wheel

Ergonomic, Attractive and Sturdy

Seat tilting system allows
the occupant to push back
and recline easily

Upper backrest with push bar
(helps the caregiver
tilt the chair)



The ideal combination of comfort and function

Providing height adjustability, versatility and excellent value, the Harmony patient care recliner can be used for taking blood samples, specific treatments, general care and prolonged sitting. It is suitable for both short and long-term care establishments.



C-C 1



C-C 2



C-C 4




C-C 6



C-C 8



C-C 9



C-C 10



C-C 11



RB-18

OPTIONS AND ACCESSORIES

ITEM	DESCRIPTION
C-C1	Left or right lateral tray
C-C2	Black adjustable elbow supports
C-C3	Adjustable IV pole, rear installation
C-C3B	Adaptor for IV pole with 1/2" round tube
C-C4	Elevating footrest support
C-C5	Set of four (4) swivel wheels with brakes
C-C6	Self-storing, removable footrests
C-C7	Side bolster cushions
C-C8	Opaque activity tray
C-C8-2	Activity tray with polycarbonate (Lexan) window
C-C9	Lateral support bracket for storing activity tray
C-C10	PRE-CURE® PRO cushion allows the insertion of speciality pressure management gels/foams
C-C11	Inclined armrest for blood sample
AB-08	Contoured headrest
AC-04	Leatherette padded tray cover
AE-02	Catheter drainage bag hook
RA-03	Lateral tray straps
RA-03A	Wrap-around tray strap
RA-19	Universal perineal abductor kit with black fixing mat
RA-23	4" (10 cm) lateral trunk supports with backrest canvas
RA-24	6" (15 cm) lateral trunk supports with backrest canvas
RB-18	Cervical support

Vinyl upholstery standard. Crypton fabric and Silvertex vinyl available as options.
Meets CA TB117-2013 fire standards.

TECHNICAL SPECIFICATIONS

C-207-20	
Overall width	28 1/2" / 72 cm
Upright length	36 1/2" / 93 cm
Reclined length	70" / 178 cm
Upright height	46 1/2" / 118 cm
Seat/floor height	Adjustable from 17 1/2" to 20 1/2" (44 to 52 cm)
Seat width	20 1/4" / 51 cm
Seat depth	Adjustable from 17" to 18 1/2" (43 to 47 cm)
Backrest height from seat	34" (86 cm)
Seat/backrest angle	Adjustable from 95 to 108 degrees
Weight w/o accessories	80 lbs. / 36 kg
Maximum weight capacity	250 lbs. / 114 kg

Please note that technical specifications are subject to change without notice.
All dimensions ± 1/2" (1.3 cm).

our complete product range

MAINTAINING PATIENT AUTONOMY



THERA-GLIDE®
WOOD



THERA-GLIDE®
METAL



THERA-GLIDE®
TRANSITIONAL



THERA-GLIDE®
OASIS ROCKING



THERA-GLIDE®
OASIS STATIONNARY



DYN-ERGO



EVOLUTION



NOMAD



VERSATILE



UTILITY

COMFORT AND POSITIONING



COMP-POSITION



HARMONY



OPTIMUM

ACCESSORIES



PRE-CURE® PRO
CUSHION



ANTI-SLIP
SURFACES

Visit our website or contact us for more details about our range of products.



lpamedical.com

2527 Dalton avenue, Quebec (Qc) G1P 3S6 CANADA
T 418-681-1313 1-800-663-4863 F 418-681-4488

Edited on 2023-01