

multi-purpose armchair

Utility

U Series

FOCUSED ON WHAT MATTERS



SPECIALTY
HEALTHCARE
SEATING



U – 120

Shown with Silvertex Imperial vinyl,
polyethylene armrest (U-C5)
and optional front wheels (U-C2)



MAINTAINING PATIENT AUTONOMY

Stackable and sturdy

Adaptable, functional and sturdy

In long-term care, armchairs of this kind are subject to intensive use. In fact, they are often moved with the patient to be closer to a table, a sink, or during an activity. To face those challenges, the UTILITY was designed to be sturdy, stackable and to be able to withstand heavy use. Whether for common areas, dining rooms or bedrooms, it will certainly meet all your needs.



TECHNICAL SPECIFICATIONS

	U-120	U-120-24	U-120-30
Overall depth	25 1/2" / 65 cm	25 1/2" / 65 cm	25 1/2" / 65 cm
Overall width	23 1/2" / 60 cm	26 1/4" / 67 cm	31" to 34" / 79 to 86 cm
Overall height	31" to 34" / 79 to 86 cm	31" to 34" / 79 to 86 cm	29 1/2" / 75 cm
Width between armrests	19 1/4" / 49 cm	22 1/4" / 57 cm	29 1/2" / 75 cm
Seat depth	17 3/4" / 45 cm	17 3/4" / 45 cm	17 3/4" / 45 cm
Seat-to-floor height	16" to 19" / 41 to 48 cm	16" to 19" / 41 to 48 cm	16" to 19" / 41 to 48 cm
Weight	25 lbs. / 11.4 kg	31 lbs. / 14 kg	N/A
Weight capacity	350 lbs. / 159 kg	450 lbs. / 204 kg	650 lbs. / 295 kg

Vinyl upholstery standard. Crypton fabric and Silvertex vinyl available as an option. Meets CA TB117-2013 fire standards.