

patient room chairs

Versatile

and **Nomad** Z Series



FOCUSED ON WHAT MATTERS

SPECIALTY
HEALTHCARE
SEATING



VERSATILE Z-25
Adjustable tilt
and backrest recline
Shown with
Z-C9 options (side panels)
and RB-18 (cervical support)



VERSATILE Z-15
Fixed seat

NOMAD Z-55

Shown with optional swing-away
footrests (option Z-C18)

MAINTAINING PATIENT AUTONOMY

Versatile and Nomad

patient room chairs

Ergonomic, comfortable, attractive and sturdy

Versatile



VERSATILE Z-25

Standard Silvertex vinyl
Available in several colours

Backrest with movable
elastic bands help balance
pressure distribution
for optimal comfort

Height-adjustable legs
(standard)



Option Z-C2-A
gas-cylinder
seat tilt

Nomade

2 directional
rear wheels



NOMAD Z-55

Fully equipped mobile
lounge chair with gas-cylinder
assisted tilt and recline


2 pivoting front wheels with brakes

A versatile and adaptable armchair

Whether stationary (Versatile) or mobile (Nomad), all standard or optional adjustments are now done easily and without tools. Its superior comfort and contemporary styling make it a patient room chair that integrates well at home or in acute or prolonged care. Its multiple sizes and options are the solution to help rehabilitation and maintain autonomy.



Z-C5



Z-C8



Z-C10



Z-C11



Z-C14



Z-C18



Z-C130R



RA-19



RA-24

OPTIONS AND ACCESSORIES

ITEM	DESCRIPTION	Widths					
		Z-15	Z-25 17	Z-25 20	Z-25 24	Z-25 28	Z-55
AE-02	Universal support for catheter drainage bag	0	0	0	0	0	0
Z-C1A	Seat height adjustment from 20 1/2" to 23 1/2" (52 to 60 cm)	0	0	0	0	0	—
Z-C2-A	Gas-spring-assisted seat tilt	—	0	0	0	+	+
Z-C3-A	Tool-less seat depth adjustment and 4-position telescopic support-operated backrest recline	—	+	0	+	+	+
Z-C4-A	Tool-less seat depth adjustment and gas-spring-assisted backrest recline	—	0	0	0	0	+
Z-C5	Rear wheels to facilitate chair transport (not to be used with a seated person)	0	0	0	0	+	—
Z-C6	Lateral support bracket for storing activity tray	0	0	0	0	0	0
Z-C6A	Rear support brackets for storing activity tray	0	0	0	0	—	0
Z-C8	PRE-CURE® PRO seat cushion	0	0	0	0	0	0
Z-C9	Upholstered side panels	0	0	0	0	0	+
Z-C10	Folding self-storing side table	0	0	0	0	0	0
Z-C11	Friction reducing leg tips	0	0	0	0	0	—
Z-C14	Opaque activity tray	0	0	0	0	0	0
Z-C15	Activity tray with Lexan window	0	0	0	0	0	0
Z-C18	Self-storing and removable footrests	—	—	—	—	—	0
Z-C130L	Left fully-adjustable elevating leg rest	0	0	0	0	0	0
Z-C130R	Right fully-adjustable elevating leg rest	0	0	0	0	0	0
RA-03	Lateral tray straps	0	0	0	0	0	0
RA-03A	Wrap-around tray strap	0	0	0	0	0	0
RA-19	Universal perineal abductor kit	0	0	0	0	0	0
RA-23	4" (10 cm) lateral trunk supports with backrest canvas	0	0	0	0	0	0
RA-24	6" (15 cm) lateral trunk supports with backrest canvas	0	0	0	0	0	0
RB-18	Cervical support	0	0	0	0	0	0

+: Standard 0: Optional —: Not available

TECHNICAL SPECIFICATIONS

	Z-15 / Z-25			Z-25-28	Z-55
	17	20	24		
Overall length	28" / 71 cm			28" / 71 cm	30 1/2" / 77 cm
Overall width	23 1/2" / 60 cm	27 1/2" / 70 cm	30 1/2" / 77 cm	33" / 84 cm	27 1/2" / 70 cm
Height	41" to 44" 104 to 112 cm	44" to 47" / 112 to 119 cm		44" to 47" 112 to 119 cm	45" / 114 cm
Seat width	17" / 43 cm	21" / 53 cm	24" / 61 cm	27 1/2" / 70 cm	21" / 53 cm
Seat depth	17 1/2" / 45 cm (Z-C3-A included)	17 1/2" / 45 cm	17 1/2" / 45 cm (Z-C3-A included)	17 1/2" / 45 cm (Z-C3-A included)	17 1/2" / 45 cm
Z-25/Z-55 seat depth (with option Z-C3-A/Z-C4-A)	Adjustable from 15 3/4" to 19 1/4" / 40 to 49 cm			Adjustable from 15 1/2" to 19 1/4" / 40 to 49 cm	Adjustable from 15 3/4" to 19 1/4" / 40 to 49 cm
Seat-to-floor height	Adjustable from 16 1/2" to 19 1/2" / 42 to 50 cm			Adjustable from 16 1/2" to 19 1/4" / 42 to 49 cm	17 1/2" / 44 cm
Weight without accessories (Z-15)	30 lbs. / 14 kg	34 lbs. / 15 kg	37 lbs. / 17 kg	N/A	N/A
Weight without accessories (Z-25/Z-55)	32 lbs. / 15 kg	36 lbs. / 16 kg	41 lbs. / 19 kg	85 lbs. / 39 kg	41 lbs. / 19 kg
Maximum weight capacity	300 lbs. / 136 kg		450 lbs. / 204 kg	750 lbs. / 340 kg	250 lbs. / 114 kg

Silvertex vinyl upholstery standard, Crypton fabric available as an option.
Meets CA-117-2013 fire standards.

our complete product range

MAINTAINING PATIENT AUTONOMY



THERA-GLIDE®
WOOD



THERA-GLIDE®
METAL



THERA-GLIDE®
TRANSITIONAL



THERA-GLIDE®
OASIS ROCKING



THERA-GLIDE®
OASIS STATIONNARY



DYN-ERGO



EVOLUTION



NOMAD



VERSATILE



UTILITY

COMFORT AND POSITIONING



COMP-POSITION



HARMONY



OPTIMUM

ACCESSORIES



PRE-CURE® PRO
CUSHION



ANTI-SLIP
SURFACES

Visit our website or contact us for more details about our range of products.



lpamedical.com

2527 Dalton avenue, Quebec (Qc) G1P 3S6 CANADA
T 418-681-1313 1-800-663-4863 F 418-681-4488

Edited on 2023-01