

Anti-slip surfaces

FOCUSED ON WHAT MATTERS

lpa+

SPECIALTY
HEALTHCARE
SEATING



ACCESSORIES

The key to good posture

Customizable, functional and durable

Extremely durable, water repellent and easy to clean: Anti-slip surfaces enable you to position a person in a chair or a bed and keep them securely in place. Our fabric strips are made of oriented fibres, so when placed one over the other, they easily slide in one direction, but not in the opposite. Our surfaces help to position the resident properly while reducing the risk of injuries to caregivers.



07-1050



07-2060



07-3050



07-2050



07-3060



07-4050

TECHNICAL SPECIFICATIONS

	07-1050	07-2050	07-2060	07-3050	07-3060	07-4050
Type	Sleeve	Sleeve	Wallet with handle	Sleeve with side handles (bariatric)	Wallet with handle	Sleeve with side handles for bed
Width	13 1/2" [34 cm]	17" [43 cm]	17" [43 cm]	21" [53 cm]	21" [53 cm]	32" [81 cm]
Length	17 1/2" [44 cm]	17 1/2" [44 cm]	52" [132 cm]	21" [53 cm]	52" [132 cm]	40" [102 cm]

Please note that these specifications are subject to change without notice.

Visit our website or contact us for more details about our range of products.



lpamedical.com

2527, avenue Dalton, Québec (Qc) G1P 3S6 CANADA
T 418 681-1313 1 800 663-4863 F 418 681-4488