

multi-purpose armchair

Utility

U Series

FOCUSED ON WHAT MATTERS



SPECIALTY
HEALTHCARE
SEATING



U – 120

Shown with Silvertex Imperial vinyl,
polyethylene armrest (U-C5)
and optional front wheels (U-C2)



MAINTAINING PATIENT AUTONOMY

Stackable and sturdy

Adaptable, functional and sturdy

In long-term care, armchairs of this kind are subject to intensive use. In fact, they are often moved with the patient to be closer to a table, a sink, or during an activity. To face those challenges, the UTILITY was designed to be sturdy, stackable and to be able to withstand heavy use. Whether for common areas, dining rooms or bedrooms, it will certainly meet all your needs.



TECHNICAL SPECIFICATIONS

	U-110	U-120	U-120-24	U-120-30
Overall depth	21 1/2" / 55 cm	25 1/2" / 65 cm	25 1/2" / 65 cm	25 1/2" / 65 cm
Overall width	22 3/4" / 58 cm	23 1/2" / 60 cm	26 1/4" / 67 cm	31" to 34" / 79 to 86 cm
Overall height	32 3/4" / 83 cm	31" to 34" / 79 to 86 cm	31" to 34" / 79 to 86 cm	29 1/2" / 75 cm
Width between armrests	18 1/2" / 47 cm	19 1/4" / 49 cm	22 1/4" / 57 cm	29 1/2" / 75 cm
Seat depth	17 3/4" / 45 cm	17 3/4" / 45 cm	17 3/4" / 45 cm	29 1/2" / 75 cm
Seat-to-floor height	18 1/4" / 46 cm	16" to 19" / 41 to 48 cm	16" to 19" / 41 to 48 cm	16" to 19" / 41 to 48 cm
Weight	19 lbs. / 8,6 kg	25 lbs. / 11,4 kg	31 lbs. / 14 kg	N/A
Weight capacity	250 lbs. / 113 kg	350 lbs. / 159 kg	450 lbs. / 204 kg	650 lbs. / 295 kg

Vinyl upholstery standard. Crypton fabric and Silvertex vinyl available as an option. Meets CA TB117-2013 fire standards.

Visit our website or contact us for more details about our range of products.



lpamedical.com

2527 Dalton avenue, Quebec (Qc) G1P 3S6 CANADA
T 418-681-1313 1-800-663-4863 F 418-681-4488